

# tac:life

respond first with confidence

the world's first  
body worn AED and  
bleed control kit



## Covert and overt options for military, medical, close protection, and civilian use.

Thanks to the **IFAB™** and **CellAED®**, it is now possible to carry an AED, as well as life-saving bleed control kit including the 'Big Four' you'll need to control traumatic haemorrhaging wherever you go.

This is unprecedented and will change the way we operate in austere and hostile environments, on the street, and in the workplace, as well as in discreet close protection settings.

The first-of-its-kind Individual First Aid Belt (**IFAB™**) gives users unparalleled peace of mind knowing that they have the core equipment required to save a life on their person, at all times.

Visit: [tac-life.co.uk](http://tac-life.co.uk) / follow: [@taclifeldtd](https://twitter.com/taclifeldtd)

**tac:life**  CorMed-DX  
respond first with confidence


Tourniquet



Chest Seal



Basic Life Support Card



Marker Pen



CellAED



Trauma Shears



Haemostatic Dressing



4" Compression Bandage



Razor (optional extra)



Our products are designed by Military Veterans and Active Close Protection professionals, and are made using high quality, Military grade materials.

CellAED® is supplied to us by our equipment partner, CorMed-DX:  
[www.cormed-dx.com](http://www.cormed-dx.com)

## The IFAB™ is made up of two main components:

### THE MAGCORE™

the mounting platform for the IFAB. It has two metal plates that line up with reciprocating magnets on the second component, the MAGCARD™. It is low profile and sits comfortably on the user's lower back/waist.

The MAGCORE™ comes in two variants:

- The MAGCORE Molle™ attaches to any existing molle pals retention system
- The MAGCORE Belt™ attaches to any belt up to two inches in width.

### THE MAGCARD™

attaches to the MAGCORE via its bespoke rare earth magnets with more than 12kg pull, and velcro. Its pouches and sleeves, made of the highest quality materials, carry the AED and all the kit you need to control the bleed.

The IFAB AED™ with full control bleed kit weighs less than a kilogram (dependant on products included), and measures 31cm by 12.7cm.





# tac:life

respond first with confidence



## F+RST AED Responder



The FirstAED responder box is designed to provide all the kit and equipment needed to tend to trauma and cardiac incidents. Complete with a bleed control kit and the revolutionary CellaAED®, all in one place, fit for purpose. The box can be mounted on walls as well as transported in vehicles or carried. This equipment is an affordable and versatile option to make sure your family remains safe. This can be wall mounted, carried, or stored in a vehicle, making sure you have the right kit there when you need it.


### The FirstAED Responder Box Includes...

- Items held in waterjet cut foam inserts
- 2 x CAT tourniquets
- 2 x Trauma Bandage 4"
- 2 x Foxseal chest seals
- 2 x Celox rapid
- 2 x pair of nitrile gloves
- 1 x CellaAED®
- 1 x Tuff cuts
- 1 x Gallant disposable razor blade
- 2 x NPA (6 & 7)
- 1 x Sharpie pen
- 1 x ATMIST/MARCH card

Graphics show how to apply the device on both adults and infants

A perforated label marks the place for easy snap activation

Audio speaker gives instructions on what to do every step of the way



# AED

## Specifications

Infant mode button for use on infants weighing less than 10kg

**Length:** 19.6cm  
**Width:** 9.3cm  
**Depth:** 1.7cm  
**Weight:** Approx 300g

LED indicator with periodic blinking showing the status of the device

# tac:life

respond first with confidence



[info@tac-life.co.uk](mailto:info@tac-life.co.uk)



0333 050 9493

SCAN TO VISIT  
[tac-life.co.uk](http://tac-life.co.uk)



Visit: [tac-life.co.uk](http://tac-life.co.uk) / follow: [@taclifeldt](https://twitter.com/taclifeldt)

**tac:life**  
respond first with confidence

 CorMed-DX